

A Novel Force Discrimination Assay Using Magnetic Beads

Chris Cole & Michael Malito

Nova Research, Inc., Alexandria, VA

Mohan Natesan*

Geo-Centers, Inc., Fort Washington, MD

*(*Protiveris, Rockville, MD)*

Rich Colton & Lloyd Whitman

Naval Research Laboratory, Washington, DC


*Supported by the Joint Science & Technology Panel for
Chemical & Biological Defense*





email: richard.colton@nrl.navy.mil

<http://stm2.nrl.navy.mil>

Quartz Sphere on AFM Cantilever



Using AFM to Measure Forces Between Single Molecules



G.U Lee, D.A. Kidwell & R.J. Colton, *Langmuir* **10**, 354 (1994).

Molecular Recognition & Force Discrimination


- Range of measured rupture force

DNA (20 mer in 100 mM NaCl)	750 +/- 120 pN
(20 mer in 10 mM NaCl)	550 +/- 70 pN
Streptavidin-Biotin	200 - 300 pN
Antigen-Antibody	40 pN

- If force discrimination is implemented in a biosensor we anticipate **increased**
 - **selectivity** (fewer false positives and negatives)
 - **sensitivity** (lower background)



Immunobead Force Discrimination Assay

- Use paramagnetic beads as labels
- Magnetic field gradient removes nonspecifically bound beads
- *Magnetic force eliminates background*




Revolutionary Biosensor Technology

Three Detection Methods




Piezoresistive
cantilever
FABS

D.R. Baselt, *et al.*, *Proc. IEEE* **85**, 672 (1997)



Transparent substrate
& optical microscope
FDB


G.U. Lee, *et al.*, *Anal. Biochem.* **287**, 261 (2000)



Magnetoresistive
elements
BARC


R.L. Edelstein, *et al.*, *Biosensors & Bioelectronics* **14**, 805 (2000)

FDB: Ovalbumin Detection on Transparent Substrates




G. Lee, *et al.*, *Analytical Biochemistry* **287**, 261 (2000)

Optical-based Force Discrimination Biosensor




magnetic bead
transparent substrate
optical microscope


Polystyrene microtiter well



Inverted optical microscope




Ultrafiltration membrane



Immunobead Force Discrimination Assay


Magnetic Beads

PEI-PEG-Antibody Chemistry



Dynal, 2.8 μm size beads containing, polymer coated mono-dispersed magnetic material (Fe_3O_4 and $\delta\text{Fe}_2\text{O}_3$)

$M \sim 12 \text{ emu/cm}^3$





SeraMag, 0.8 μm size beads containing a magnetite core coated with polymer

$M \sim 38 \text{ emu/cm}^3$

Immunobead Force Discrimination Assay

Results for Ovalbumin



- Detection sensitivity = 10 pg/mL
- Assay time = 25 min.

Immunobead Force Discrimination Assay

Conclusions

- Sensitivity of our immunobead force discrimination assay is $\sim 10^3$ X higher than conventional immunoassays such as ELISA
- Specificity is typically $> 99\%$
- Response time is comparable to PCR-based instruments
- Better magnetic beads are needed to improve force discrimination between molecules